# **Retinopathy Detection**

AI client project 3

**Brief Documentation** 

#### 1. Introduction

This project involves building a machine learning model to detect the type of diabetic retinopathy from retinal images. The final model is deployed as a web application using Flask and hosted on an open-source platform (Render).

### 2. Dataset

The dataset contains labeled retinal images classified into three categories:

- No DR (No Diabetic Retinopathy)
- Mild DR
- Severe DR

The dataset was used to train, validate, and test the deep learning model.

## 3. Data Preprocessing

Images were resized to 224x224 pixels, normalized, and augmented. The images were converted into NumPy arrays and split into training and validation datasets.

## 4. Model Building

A Convolutional Neural Network (CNN) with a custom focal loss function was used. The model was compiled with Adam optimizer and categorical crossentropy/focal loss as the loss function. The best performing model was saved using model checkpointing.

#### 5. Model Evaluation

The model was evaluated using metrics such as accuracy, precision, recall, and F1-score. The classification report showed high performance, especially for the 'No DR' and 'Severe DR' classes.

## 6. Web Application

A Flask-based web app was built where users can upload an image and get the predicted retinopathy type. The app loads the trained model and returns predictions after processing the image.

## 7. UI Design

The HTML template was styled to provide a professional and clean interface. Users can upload images and view the result instantly.

## 8. Deployment on Render

The project was deployed using Render. The steps included:

- Creating a GitHub repository and pushing all files.
- Adding `requirements.txt` and `render.yaml` files.
- Linking the GitHub repo with Render and deploying the app.
- Testing the deployment and ensuring predictions work on uploaded images.

#### 9. Conclusion

This project demonstrates how deep learning can be effectively used for medical image classification. It also highlights how to deploy a model as a web application on a cloud platform. You can access it on https://retinopathy-g5tz.onrender.com